







REF 46315 	ZENIT RA LKM-1	Distributed by 
INSTRUCTIONS FOR USE	   50	

USAGE

The *ZENIT RA LKM-1 (Liver Kidney Microsomes type 1)* test is a chemiluminescent immunoassay (CLIA) used for quantitative determination, using dedicated *ZENIT RA Analyser* instrumentation, of specific IgG class antibodies directed against Liver Kidney Microsomes type 1 in samples of human serum or plasma (EDTA, Sodium Citrate).

This dosage is used as a diagnostic aid when assessing type 2 Autoimmune Hepatitis (EAI or AIH) infections.

CAUTION: Medical decisions must not be based exclusively on the result of this test, but must take into account all available clinical and laboratory data as a whole.

CLINICAL SIGNIFICANCE

Autoimmune Hepatitis (AIH) is a chronic inflammation of the liver of unknown origin, in which immune tolerance to hepatocytes is lost¹. The disease is detected by the presence of a number of characteristic alterations from the histological point of view, when blood tests show the presence of gamma globulins and antibodies². The disease is found mainly in females, both of paediatric and adult ages³. It is an autoimmune process and therefore is associated, sometimes, with other autoimmune disorders such as ulcerative colitis⁴ and Graves' Disease⁵.

Autoimmune Hepatitis is classified according to the specific autoantibodies detected, which are, for **type 1**, principally represented by ANA (anti-nuclear) and SMA (anti-smooth muscle) antibodies; **type 2** on the contrary is typified by LKM-1 (type 1 liver kidney microsomal) antibodies; **type 3** is characterised by anti-SLA (soluble liver antigen) and anti-LP (liver pancreas antigen) antibodies.

The presence of anti-LKM autoantibodies was described for the first time in 1973⁶. It was subsequently demonstrated that these are a heterogeneous family of autoantibodies and, depending on the target antigen, are classified in three types known as anti-LKM 1, 2 and 3. Cytochrome P4502D6 was identified as the autoantigen of the anti-LKM-1 antibodies⁷. Cytochrome P4502C9 and UDP-glucuronosyltransferase were later identified respectively as the targets of the anti-LKM-2 and anti-LKM-3 antibodies⁸. While anti-LKM-2 antibodies have been associated with drug-induced hepatitis, anti-LKM-1 antibodies have been associated with **type 2** AIH and are a characteristic marker for this disease, not usually associated with the presence of ANA and/or ASMA of **type 1** AIH.

However, anti-LKM-1 antibodies are not specific to **type 2** AIH, in that they may be found in 5-10% of the patients with chronic HCV hepatitis⁹, even though the epitopes recognised in HCV patients are different from those recognised in patients suffering from **type 2** AIH^{10,11}.

Anti-LKM-1 antibodies are found in around 90% of patients with **type 2** AIH⁸ and can be identified by the IFA test using sections of mouse liver and kidney tissue. It is also possible to test for anti-LKM-1 antibodies with indirect solid phase assays, using the recombinant protein CYP2D6 with high levels of sensitivity and specificity^{12,13}.

PRINCIPLE OF THE METHOD

The *ZENIT RA LKM-1* kit for quantitative determination of specific IgG class antibodies directed against Liver Kidney Microsomes employs an indirect two-step immunological method based on the principle of chemiluminescence.

The specific antigen used to coat the magnetic particles (solid phase) and the G protein able to bond with human IgGs is marked with an acridinium ester derivative (conjugate).

During initial incubation, the specific antibodies present in the sample, in the calibrators or in the controls bond with the solid phase.

During the second incubation, the conjugate reacts with the anti-LKM-1 antibodies sequestered by the solid phase.

After each incubation, the material that has not bonded with the solid phase is removed by aspiration and subsequent washing.

The quantity of marked conjugate that remains bonded to the solid phase is assessed by activation of the chemiluminescence reaction and measurement of the luminous signal. The generated signal, expressed in relative light units (RLU), is indicative of the concentration of specific antibodies present in the sample, in the calibrators and in the controls.

AUTOMATION

The *ZENIT RA Analyser* instrument automatically performs all the operations envisaged by the dosage protocol: addition of samples, calibrators, controls, magnetic particles, conjugates and chemiluminescence activation solutions to the reaction container; magnetic separation and washing of particles; measurement of the light emitted.

The system calculates the dosage results for the samples and controls by means of a stored calibration curve and prints a report that includes all the information related to the dosage and to the patient.

MATERIALS AND REAGENTS

Materials and reagents supplied

REAG	1	MP	2.5 mL
------	---	----	--------

Magnetic particles coated with LKM 1/cytochrome P450 2D6 recombinant antigen in Phosphate Buffer containing stabilising proteins, surfactant, Pro-Clin 300 and sodium azide (< 0.1 %) as preservatives.

REAG	2	DILCONJ	10 mL
------	---	---------	-------

Conjugated Diluent Solution: Phosphate Buffer containing stabilising proteins, surfactant, Pro-Clin 300 and sodium azide (< 0.1 %) as preservative.

REAG	3	CONJ	4 mL
------	---	------	------

G protein conjugated with an acridinium ester derivative (conjugate), in Phosphate Buffer containing stabilising proteins, surfactant, Pro-Clin 300 and sodium azide (< 0.1%) as preservative.

REAG	4	DIL	25 mL
------	---	-----	-------

Sample Diluent Solution: Phosphate Buffer containing bovine serum albumin, a surfactant, an inert blue colouring agent, Pro-Clin 300 and Gentamicin SO₄ as preservatives.

REAG	5	CAL A	0.5 mL
------	---	-------	--------

Serum with low concentration of anti LKM-1 antibodies containing sodium azide (< 0.1%) as preservative.

REAG	6	CAL B	0.5 mL
------	---	-------	--------

Serum with high concentration of anti LKM-1 antibodies containing sodium azide (< 0.1%) as preservative.

All reagents are ready for use.

Reagents 1, 2, 3 and 4 are assembled in a single kit forming the reagent cartridge.

The Calibrator concentrations are expressed in AU/mL (arbitrary units) and calibrated against an internal reference standard. The concentration settings, specific for each production batch, are recorded on the DATA DISK included in the kit.

DATA DISK

A Mini-CD containing data regarding all the products in the ZENIT RA line (Reagents, Calibrators, Control Serums) updated to the last production batch with the exclusion of products that have expired at the date when the new DATA DISK was compiled.

Only the DATA DISK with the highest batch number needs to be kept to maintain the information required for correct operation of the system up to date.

Materials and reagents required but not supplied in the kit:

- | | |
|---|----------------|
| - ZENIT RA Analyser ⁽¹⁾ | Code No. 41400 |
| - ZENIT RA Cuvette Cube ⁽¹⁾
Pack of 960 cuvettes | Code No. 41402 |
| - ZENIT RA System Liquid ⁽²⁾
1 bottle containing 5 litres of ready-to-use solution | Code No. 41409 |
| - ZENIT RA Wash Solution ⁽²⁾
1 bottle containing 10 litres of ready-to-use solution | Code No. 41407 |
| - ZENIT RA Trigger Set ⁽²⁾
1 250 mL-bottle of Trigger A (pre-trigger solution)
1 250 mL-bottle of Trigger B (trigger solution) | Code No. 41403 |
| - ZENIT RA D-SORB Solution
Pack of 2 bottles containing 1 litre of ready-to-use solution | Code No. 41436 |
| - ZENIT RA Cartridge Checking System ⁽²⁾ | Code No. 41401 |
| - ZENIT RA Top Cap Set
300 red top caps to close the calibrator containers after first use | Code No. 41566 |

⁽¹⁾ Made by IDS France SAS, 42 rue Stéphane Mazeau, 21320 Pouilly en Auxois, France and distributed by A. Menarini Diagnostics Srl.

⁽²⁾ Made by IDS SA, 101-103 rue Ernest Solvay, 4000 Liège, BELGIUM and distributed by A. Menarini Diagnostics Srl.

Other Recommended Reagents

ZENIT RA LIVER CONTROL SET

Code No. 46317

Three 1.0 mL vials of negative human serum and three 1.0 mL vials of human serum positive for anti-LKM-1 antibodies.

WARNINGS AND PRECAUTIONS

The reagents supplied in the *ZENIT RA LKM-1* kit are only for in vitro diagnostic use and not for in vivo use in humans or animals.

This product must be used in strict compliance with the instructions given in this document by professional users.

Menarini cannot be held responsible for any losses or damages caused by use not in conformity with the instructions supplied.

Safety precautions

This product contains material of animal origin and therefore must be handled as if it contains infecting agents.

This product contains components of human origin. All units of serum or plasma used to produce the reagents in this kit have been analysed with FDA-approved methods and found not to be reactive due to presence of HBsAg, anti-HCV, anti-HIV1 and anti-HIV2.

However, since no analysis method is able to guarantee the absence of pathogenic agents, all material of human origin must be considered to be potentially infected and handled as such.

In the event of damaged packaging or accidental leakage, decontaminate the area concerned with a diluted solution of sodium hypochlorite after putting on suitable personal protective equipment (overall, gloves, goggles).

Dispose of the material use for the clean-up and of the packaging involved in the leakage according to national regulations for disposal of potentially infected waste.

Some reagents contain sodium azide as a preservative. Since sodium azide may react with lead, copper and leaded brass forming explosive azides in piping, it is recommended that reagents or waste are not poured down drains but are disposed of in compliance with the national regulations on disposal of potentially hazardous waste.

Operating precautions

Reliable results can only be obtained by strictly complying with these instructions and scrupulously following what is written in the operating manual for the instrument.

The reagents supplied in the kit must be used only with the ZENIT RA *Analyser* system.

The components of the reagent cartridge must not be removed from the cartridge and reassembled.

Do not use the kit after its use-by date.

PREPARATION OF THE REAGENTS

The reagents supplied in the kit are all ready for use.

PRESERVATION AND STABILITY OF THE REAGENTS

Store the reagents supplied in the kit at 2-8°C in a vertical position in a dark place.

In these conditions unopened reagent cartridge and calibrators are stable up to the use-by date.

After opening the reagent cartridge can be used for 60 days if kept in a refrigerator at 2-8°C or in the analyser.

After opening, the calibrators can be used for 60 days if kept in a refrigerator at 2-8°C and if they have not been left in the analyser for more than 6 hours per session.

Do not freeze the reagents and calibrators.

PREPARATION AND PRESERVATION OF SAMPLES

Dosage must be performed on samples of human serum and plasma (EDTA – Sodium Citrate).

Use of lipaemic, haemolysed and turbid samples is not recommended.

If the dosage is performed after more than 8 hours, separate the serum or the plasma from the clot, from the red globules and from the separating tubes with gel.

Before being analysed, samples may be kept in a refrigerator at 2-8°C for a maximum of 7 days.

If the dosage is to be performed after more than 7 days, store the samples frozen (< - 20°C).

Avoid repeated freezing and thawing.

OPERATING PROCEDURE

Scrupulously follow the instructions given in the operating manual of the instrument to obtain reliable analytical results.

Loading of reagents

All the reagents supplied in the kit are all ready for use.

Before inserting the reagent cartridge in the system, the magnetic particle container must be horizontally agitated by rotation in order to favour resuspension of the particles. Avoid generating foam when performing this operation.

Place the reagent cartridge in the reagent area of the instrument using the guide provided and leave it to be agitated for at least 40 minutes before use.

Positioning of the reagent cartridge simultaneously determines reading of the identity bar-code. If the cartridge label is damaged or if it is not readable, the reagent cartridge identification data can be entered manually.

The instrument automatically maintains the magnetic particles constantly agitated.

If the reagent cartridge is removed from the instrument, store it at 2-8°C in a vertical position in a dark place.

Loading of calibrators

ZENIT RA calibrators are ready for use. Leave the calibrators at room temperature for 10 minutes and then gently shake the contents, either manually or using a vortex, avoiding the formation of foam.

When using the calibrators for the first time, remove the guarantee seal and the white sealing cap before placing them in the analyser.

If the calibrators have already been used, the container will have a top cap (red cap) with no guarantee seal. Remove the red closing cap before placing them in the analyser.

Place the calibrators in the samples area of the analyser; see the analyser user manual on how to identify them in the analyser. Bar-code data must be entered manually if the label is damaged or if it is unreadable.

The readings for the LKM-1 antibody concentration in the calibrators are recorded in the DATA DISK and automatically transferred to the analyser.

At the end of the session, the calibrator containers must be closed with the top caps (red caps) provided and stored at 2-8°C until they are used again.

The calibrators can be used for a maximum of four times.

Loading of controls

Place the controls in the samples area of the analyser. See the analyser user manual on how to identify them in the analyser. If there is no bar-code on the control or if it is not readable, the control identification data must be entered manually. If Zenit RA Controls are used, see the usage instructions provided. The readings for the LKM-1 antibody concentration in the Zenit RA controls are recorded in the DATA DISK and automatically transferred to the analyser.

Select the required parameters for each control.

Loading of samples

Place the samples in the samples area of the analyser; see the analyser user manual on how to identify them in the analyser. If there is no bar-code on the sample or if it is not readable, the sample identification data must be entered manually.

Select the required parameters for each sample.

Calibration

The *ZENIT RA Analyser* instrument uses a memorised calibration curve (master curve), generated by the manufacturer for each batch of reagent cartridges.

The “master curve” parameters, together with the calibrator concentration settings, are stored in the DATA DISK and transferred to the instrument’s data base.

Calibrators A and B are used to recalibrate the “master curve” in both for the instrument used and for the reagents on board.

To perform recalibration analyse the two calibrators A and B in triplicate and the controls singly. The concentration readings obtained with the controls make it possible to validate the new calibration.

Once recalibration of the “master curve” has been accepted and memorised, all subsequent samples can be analysed without any further calibration, except in the following cases:

-
- when a reagent cartridge with a new batch number is loaded into the instrument;
 - when the control readings do not fall within the range of acceptability;
 - when the instrument maintenance procedure is performed;
 - when the validity of the recalibrated “master curve” has expired.

The validity of the recalibrated “master curve” for the *ZENIT RA LKM-1* kit is 21 days.

Control of recalibration is performed automatically by the instrument.

Dosage

Press the start button.

1. The system aspirates 80µL of Sample Diluent, 40µL of Magnetic Particles, 100µL of Sample Diluent and 10µL of the sample, calibrators or control. The aspirated solutions and suspension are dispensed into the reaction cuvette.
2. The reaction cuvette is incubated in the rotor at 37°C for 10 minutes.
3. After this phase of incubation, the magnetic particles are separated and washed.
4. 40µL of conjugate and 160 µL of conjugate diluent are dispensed into the cuvette.
5. The reaction cuvette is incubated in the rotor at 37°C for 20 minutes.
6. After this last phase of incubation, the magnetic particles are separated and washed and the cuvette is transferred to the reading chamber.
7. The quantity of conjugate bonded to the solid phase, expressed in RLU, is directly proportionate the concentration of LKM-1 antibodies present in the sample.
8. The responses obtained are interpolated on the calibration curve and transformed into concentrations.

QUALITY CONTROL

To ensure the validity of the dosage, control serums at differing levels of concentration (at least one negative serum and one positive serum) must be measured every day in which dosage is performed.

If your laboratory requires a more frequent use or a higher number of controls to check the dosage results, comply with the set quality control procedure.

If ZENIT RA control serums are used, the expected average readings and the acceptability limits are those given on the DATA DISK included in the control pack too.

If different control serums are used, before using them, the readings expected with ZENIT RA reagents and system must be defined.

If the control reading does not fall within the specified range of acceptability, the related dosage results are not valid and the respective samples must be analysed again.

In this case, before repeating the dosage, a recalibration procedure must be performed.

CALCULATION AND INTERPRETATION OF THE RESULTS

Calculation of the results

The concentration of the anti-LKM-1 antibodies present in the samples that are being tested is automatically calculated by the system. The readings can be viewed on the display or printed.

The concentrations are expressed in AU/mL.

Calculation of the analyte concentration in the sample takes place by reading the response obtained for each sample on a calibration curve processed by a logistic “fitting” system with four parameters (4PL, Y weighted), periodically corrected according to the responses obtained for dosage of the calibrators.

For detailed information on how the system calculates the results, please see the operating manual for the system.

Interpretation of the results

The dosage range measurable for *ZENIT RA LKM-1* is: 0.0 – 300 AU/mL.

Readings lower than 0.0 AU/mL are extrapolated values, the message “OMR-” and/or ORA appears and they are shown as “equal to 0.0 AU/mL”.

Readings higher than 300 AU/mL are accompanied by the message “OMR+” and/or ORA and may be retested after suitable dilution.

The results of the samples may be interpreted in the following way:

(AU/mL)	Interpretation
< 15	The sample must be considered to be negative for the presence of anti-LKM-1 antibodies
≥ 15	The sample must be considered to be positive for the presence of anti-LKM-1 antibodies

The above readings are to be considered only suggested readings. Each laboratory must establish its own reference ranges.

DOSAGE LIMITS

For diagnostic purposes, the results obtained with the *ZENIT RA LKM-1* kit and the *ZENIT RA Analyser* system must be used together with the other clinical and laboratory data available to the physician.

Bacterial contamination of the sample and heat inactivation may influence the result of the dosage.

Heterophyllous antibodies present in human serum samples may react with immunoglobulin-based reagents, causing interference with in vitro immunological dosages. Such samples may give rise to anomalous readings if analysed with the *ZENIT RA LKM-1 kit*.

EXPECTED READINGS

The samples of 100 individuals selected randomly from the normal routine work of the laboratory were analysed to check the presence of anti-LKM-1 antibodies.

All samples analysed proved negative, with an average reading of 0.9 AU/mL and a standard deviation of 5.6 AU/mL.

Using the results obtained the "Limit of Blank" (LoB = the highest reading that can be expected in a series of samples that do not contain the analyte) was calculated. The "Limit of Blank", determined as 95° per centile of the negative population, proved equal to 1.6 AU/mL with the Reagent batch no. 2.

DIAGNOSTIC SENSITIVITY AND SPECIFICITY (CLINICAL)

A total of 423 samples were tested with the ZENIT RA *LKM-1* kit: 8 of the samples were from patients suffering from type 2 autoimmune hepatitis (AIH-2), 58 from patients suffering from type 1 autoimmune hepatitis (AIH-1), 118 from patients with primary biliary cirrhosis (PBC), 40 from patients with hepatitis C (HCV), 19 from patients with hepatitis B (HBV), 20 samples from patients with primary sclerosing cholangitis (PSC), while 160 samples were presumed normal, originating from routine laboratory testing.

In the presumably negative population studied, including 58 samples from type 1 autoimmune hepatitis (AIH-1), 118 from patients with primary biliary cirrhosis (PBC), 40 with hepatitis C (HCV), 19 with hepatitis B (HBV), 20 primary sclerosing cholangitis (PSC) and 160 normal samples, 7 samples tested positive and 408 negative.

- **Diagnostic specificity: 98.3% (408/415)**

Of the 7 samples that proved to be "not negative", 4 belonged to the group of patients diagnosed with AIH-1, 1 belonged to the group of patients diagnosed with HCV, 1 belonged to patients diagnosed with HBV, and 1 to the group of normal samples.

In the presumably positive population studied, including 8 samples from patients diagnosed with type 2 autoimmune hepatitis (AIH-2), 8 proved to be positive.

- **Diagnostic sensitivity: 100.0 % (8/8)**

Based on the diagnostic specificity and sensitivity results, **the diagnostic agreement is 98.3 % (416/423)**.

PERFORMANCES

Caution: the data presented do not represent the operating specifications of the kit, but serve as experimental proof of how the kit works within these specifications in the manner envisaged by the manufacturer.

Precision and Reproducibility

The precision and the reproducibility of the ZENIT RA *LKM-1* kit have been assessed using a protocol based on the guidelines given in Clinical and Laboratory Standards (CLSI) document EP5-A2.

The **precision** was calculated by analysing the results of 20 replicates of four serums (one negative and three positive with differing concentrations of LKM-1 antibodies) performed with two different batches of reagents in the same test session.

The concentration of the negative anti-LKM-1 serum (N1) was always 0.0 AU/mL with reagent batch no. 1 and within the range from 0 to 0.7 AU/mL with batch no. 2.

The table shows the results obtained with the 3 positive serums.

Sample	Reagents Batch no.	Average Concentration (AU/mL)	SD	CV %
1	1	65.6	2.65	4.0
	2	61.8	2.47	4.0
2	1	163.0	8.88	5.4
	2	162.0	6.00	3.7
3	1	210.7	10.73	5.1
	2	198.9	11.51	5.8

The **reproducibility** was calculated by analysing the results of the determination of five serums (one negative and four positive with differing concentrations of anti-LKM-1) performed singly, in 30 different sessions, with two different batches of reagents.

The concentration of the negative anti-LKM-1 serum (N1) was always within the range from 0.0 and 0.9 AU/mL.

The table shows the results obtained with the 4 positive serums.

Sample	Average Concentration (AU/mL)	SD	CV %
1	58.9	3.16	5.4
2	149.2	8.89	6.0
3	192.8	13.62	7.1
4	86.8	4.47	5.1

Analytical Sensitivity

The analytical sensitivity of the *ZENIT RA LKM-1* kit was assessed using a protocol based on the guidelines given in Clinical and Laboratory Standards (CLSI) document EP17-A.

In one case referred to as the **Limit of Detection** (*LoD*: that is the smallest quantity of analyte that the method is able to measure) the formula for calculating $LoD = LoB + C_{\beta} SD_s$ (in which *LoB* is the "Limit of Blank", SD_s is the estimated standard deviation of the distribution of the sample at low concentration and C_{β} is derived from 95^o percentile of the standard Gaussian distribution) was applied.

Three low concentration samples of analyte were used, determined singly with one batch of reagents in 15 different tests.

The Limit of Detection of the *ZENIT RA LKM-1* kit proved to be 3.6 AU/mL.

In the other protocol, calculation of the **Minimum Detectable Concentration** (*MDC*) is envisaged: 20 replicates of the solution of the 0 AU/mL Standard of the Master curve were used. The average and the

standard deviation (DS) were formulation on a batch of the kit. The RLU value related to the average + 2.6 DS was interpolated on the curve and the related concentration was obtained from this value.

The analytical sensitivity expressed as the Minimum Detectable Concentration is 0.7 AU/mL.

The minimum detection values, together with considerations of a clinical kind and the results of comparison with reference methods contributed to the definition of the cut-off value.

Analytical Specificity: Interferences

A study based on the guidelines given in the CLSI document EP7-A2 has shown that the dosage performances are not influenced by the presence in the sample of the potentially interfering substances listed in the table below, up to the tested concentration.

Potentially Interfering Substances	Maximum tested concentration
Free bilirubin	20 mg/dL
Conjugated bilirubin	20 mg/dL
Haemoglobin	1000 mg/dL
Fatty acids	3000 mg/dL

Use of lipaemic, haemolysed and turbid samples is not in any case recommended.

Analytical Specificity: Cross reactions

In order to assess potential cross reactions of the antigen used to sensitise the microparticles, a study was conducted using 40 samples, all with high levels of other antibodies and negative for anti-LKM IgG.

The samples used were subdivided as follows: SS-A (2), SS-B (2), U1-snRNP (2), Jo-1 (2), Scl-70 (2), Cenp B (2), Sm (2), PR3 (2), MPO (2), β_2 -GLI/CL (2), Gliadin (2), t-TG (2), CCP (2), GBM (2), dsDNA (2), TG (3), TPO (3) and Rheumatoid Factor (RF) (4).

The study did not show any significant cross reaction to the antigen in the solid phase with the other antibodies.

Saturation effect at high doses

Some immunological methods used to determine samples containing the analyte at extremely high concentrations may supply apparent levels of underestimated analyte (Hook effect).

The method used in the *ZENIT RA LKM-1* kit, since it uses two incubations, is not subject to this effect.

A sample with an extremely high concentration (above the measurement range) of anti-LKM-1 antibodies confirmed the absence of the "hook" effect up to the concentration of 22000 AU/mL.

Relative Sensitivity and Specificity

The presence of anti-LKM-1 antibodies was determined using the *ZENIT RA LKM-1* kit and an ELISA dosage method available on the market in 343 samples: 8 of the samples were from patients suffering from

type 2 autoimmune hepatitis (AIH-2), 118 from patients suffering from primary biliary cirrhosis (PBC), 58 from patients with type 1 autoimmune hepatitis (AIH-1), 40 from patients with hepatitis C (HCV), 19 from patients with hepatitis B (HBV), 20 samples from patients with primary sclerosing cholangitis (PSC), while 80 presumably normal samples were obtained from routine laboratory testing.

8 samples gave rise to discordant results between the ZENIT RA dosage and the ELISA dosage available on the market.

The **relative concordance** was therefore found to be 97.7 % (95% Confidence Interval: 95.3–98.9% (335/343)).

The **relative sensitivity** was therefore found to be 100.0 % (95% Confidence Interval: 51.7–98.5% (6/6)).

The **relative specificity** was therefore found to be 97.6% (95% Confidence Interval: 95.2–98.9% (329/337)).

BIBLIOGRAPHY

1. Strassburg CP. Autoimmune hepatitis. *Best Pract Res Clin Gastroenterol*, 2010. 24(5): 667-82.
2. Kachar R. and Fallon M. Diagnostic criteria for autoimmune hepatitis: what is the gold standard? *Hepatology*, 2010. 51(1): 350-1
3. Ngu JH et al. Population-based epidemiology study of autoimmune disease: a disease of older women? *J Gastroenterol Hepatol*, 2010. 25(10): 1681-6.
4. Navaneethan U. and Shen B. Hepatopancreatobiliary manifestations and complications associated with inflammatory bowel disease. *Inflamm Bowel Dis*, 2010. 16(9): 1598-619.
5. Nakamura H. et al. Prevalence of interrelated autoantibodies in thyroid diseases and autoimmune disorders. *J Endocrinol Invest*, 2008. 31(10): 861-5.
6. Rizzetto M, Swana G, Doniach D. Microsomal antibodies in active chronic hepatitis and other disorders. *Clin Exp Immunol* 1973; 15: 331-44.
7. Gueguen M, Meunier-Rotival M, Bernard O, Alvarez F. Anti-liver kidney microsome antibody recognizes a cytochrome P450 from IID subfamily. *J Exp Med* 1988; 168: 801-6.
8. Philipp T, Durazzo M, Trautwein C, Alex B, Straub P, Lamb JG, et al. Recognition of uridine diphosphate glucuronisyl transferases by LKM-3 antibodies in chronic hepatitis D. *Lancet* 1994; 344:578-81.
9. Lenzi M, Ballardini G, Fusconi M, Cassani F, Selleri L, Volta U, et al. Type 2 autoimmune hepatitis and hepatitis C virus infection. *Lancet* 1990; 335:258-9.
10. Declos-Valee JC, Hajaoui O, Yamamoto AM, Jacqz-Aigrin E and Alvarez F. Conformational epitopes on CYP2D6 are recognized by Liver/kidney microsomal antibodies. *Gastroenterology* 108: 467-470.1995.
11. Dalekos GN, Wedemeyer H, Obermayer-Straub P, Kayser A, Barut A, Frank H and Manns MP. Epitope mapping of cytochrome P450 2D6 autoantigen in patients with chronic hepatitis C during alpha-interferon treatment. *J. Hepatol.* 30(3): 366-375.1999.
12. Miyakawa H, Kikazawa E, Abe K, Kikuchi K, Fujikawa H, Matsushita M, Kawaguchi N, Morizane T, Ohya K, Kako M. Detection of anti-LKM-1 (anti-CYP2D6) by enzyme-linked immunosorbent assay in adult patients with chronic liver diseases. *Autoimmunity*, 1999;30(2):107-14

-
13. Kekar N, Ma Y, Davies ET, Cheeseman P, Mieli-Vergani G, Vergani D. Detection of liver kidney microsomal type 1 antibody using molecularly based immunoassay. *Journal of Clinical Pathology* 2002;55(12): 906-909.

**TECHNOGENETICS S.r.l.**

Via Vanvitelli, 4
20129 – Milan – Italy

Operating facilities: Via della Filanda, 26
26900 – Lodi - Italy

GREAT BRITAIN**Distributed by**

A. Menarini Diagnostics Ltd.
405 Wharfedale Road
Winnersh – Wokingham – Berkshire RG41 5RA
Tel. +44 118 94 44 100 - Fax +44 118 94 44 111
www.menarindiag.co.uk

NETHERLANDS**Distributed by**

Menarini Diagnostics Benelux N.V.
De Haak, 8 - 5555 XK Valkenswaard
Tel. +31 40 20 82 000 - Fax +31 40 20 42 184
www.menarindiagnosics.nl